

ABOUT THE ANTIPARASITIC ZAPPER (PRIOR CHECKS)

Contraindications:

- ✓ Pregnant or breastfeeding women;
- ✓ Pacemaker or defibrillator;
- ✓ Chemotherapy or radiation in progress;
- ✓ Dialysis/Severe renal impairment;
- ✓ Medicated chronic epilepsy;
- ✓ **Heart problems;**
- ✓ **Metal prostheses.**

Important notes:


- ✓ Hydrate well before, during, and after each session (2 liters of healthy water/day);
- ✓ Avoid heavy and high-calorie foods throughout the protocol;
- ✓ Suspend the session in case of intense and persistent discomfort;
- ✓ Avoid touching the copper rods with open and bloody wounds;
- ✓ Avoid late sessions if it gives you a strong boost of vitality;
- ✓ Clean the main unit with a damp cloth without any chemicals;
- ✓ Never submerge the main unit in water and handle it gently;
- ✓ **For certain situations**, such as young children under 6 years old, the elderly, cases of frail health, or multiple medications, avoid discomfort by starting more gently with a few sessions of only 5 minutes;
- ✓ **Vary the frequency from session to session** to avoid habituation and obtain a broader spectrum of effects on the countless different parasites;
- ✓ **Never use the Zapper while your feet are in the water** with the ionizer and the bracelet on you. Always remove accessories from the ion bath before using the Zapper;
- ✓ **The heating belt is optional** and can be used independently or simultaneously with the other functions of the 3-in-1 VitalBody;
- ✓ **The ionic bath and the Zapper do not work simultaneously;**
- ✓ The Zapper can be very useful for pets.

HOW TO USE ZAPPER?

A. COPPER RODS SET UP



- ✓ **Strongly pull the 2 ends of the main wire, slowly tearing its sheath** as much as necessary to distance the red connector from the black connector;
- ✓ Separate the wires in the opposite direction and hold them firmly under the 2 connectors;
- ✓ Connect the 2 twin ends of the wire to the ends of each tube without taking into account the red or black color of these;
- ✓ Connect the single end of the wire of the copper rods in front of the main unit;
- ✓ **Wet the skin and the copper rods** to obtain better electrical conduction;
- ✓ Install the rods knowing that the frequencies pass from one rod to the other by the shortest path, so one rod emits the signal and the other receives it;
- ✓ Position the rods opposite each other in the hands, under the feet, on the hips, under the armpits, etc.



B. STAND-ALONE VITALZAPPER OPERATING



- ✓ Read the contraindications and important notes beforehand in order to respect them;
- ✓ Connect and switch on the device using the small electrical conversion block before activating the "ON/OFF" switch;
- ✓ **The internal 9 V battery takes over** if the device is disconnected from the mains;
- ✓ Select by connection the desired frequency on the main device:
 - 30 kHz = 30,000 Hz = Dr. Hulda's High Frequency (no feeling),
 - 2.5 kHz = 2,500 Hz = RIFE intermediate frequency (few sensations),
 - 15 Hz = Low Croft frequency (perceptible sensations).
- ✓ **Time the duration of the session to 10 minutes**, which is sufficient in most cases;
- ✓ If everything is in order, a faint constant or flickering white light appears on the top of the device aligned with the connection of the wires of the copper rods;
- ✓ **In the event of a problem, follow the troubleshooting procedure below.**

C. STAND-ALONE VITALZAPPER TROUBLESHOOTING

If the device is on, but nothing happens during the session:

- ✓ Ensure that the skin in contact with the copper remains wet;
- ✓ It is normal not to feel the passage of electric current when using the 30 kHz (high) frequency and even that of 2.5 kHz (intermediate);
- ✓ Some users do not feel any of the three frequencies, such as electricians and people on antidepressants, which prevent a nervous response;
- ✓ Select the frequency of 15 Hz (low) to better feel the current's passage.

If the main unit appears to have no power:

- ✓ Tighten the electrical conversion block connections;
- ✓ Check that the **green operating light** in the upper right corner is on;
- ✓ Change the 9 V battery in the compartment screwed under the device when it is used in autonomous mode, i.e. disconnected from the electrical network.

If a problem persists, please contact us at 1-866-981-0357.

D. END OF SESSION

- ✓ **Turn off the device and unplug the power supply** from the wall outlet to avoid it being unnecessarily powered for several days;
- ✓ Rinse, wash and dry the skin that has been in contact with copper;
- ✓ If possible, leave the connectors in place at the end of each rod;
- ✓ **Drink a large glass of water with 20 drops of ConcenTrace mineral supplement** and add a little lemon or orange juice to hide the taste if necessary.

E. COPPER RODS MAINTENANCE

- ✓ It is normal and not harmful for copper to oxidize with skin acidity and become slightly brown or greenish with time and use;
- ✓ Ideally, do nothing and let the copper oxidize, which affects only its aesthetic appearance, not its effectiveness;
- ✓ Clean the rods with a magic sponge or a mild metal cleaner if desired.

F. ZAPPER USAGE PROTOCOL

- ✓ **Perform a dozen sessions with the ionizing bath beforehand for the 1st month** before using the Zapper to cleanse the body while promoting better circulation of fluids and greater filter organs' efficiency;
- ✓ Skip the previous step if you do not have the 3-in-1 VitalBody device or contact us for more information about it;
- ✓ **Run 10 minutes of Zapper daily** in the 2nd month of the protocol;
- ✓ **Run 3 x 10 minutes of Zapper daily** in the 3rd month of the protocol;
- ✓ **The Zapper is used daily as an antibiotic;**
- ✓ Start with a few 5-minute sessions only for a case of frail health risking initial detoxifying discomfort;
- ✓ **Regularly vary frequency (15 Hz, 2.5 kHz, 30 kHz) and rods position** to avoid habituation and to obtain a broader spectrum of effectiveness on parasites;
- ✓ **Follow the protocol well to avoid a toxic overload** recirculating in the lymphatic system acting as a collector of cellular waste;
- ✓ **In the event of persistent discomfort or recovery crisis** lasting more than three days, stop the sessions for a certain time to rest the metabolism;
- ✓ **Drink plenty of pure water, eat healthily, and exercise more** to promote tissue and body fluids detoxifying;
- ✓ **Resume the protocol where it was interrupted** as soon as the temporary complications and the symptoms of metabolic cleansing have disappeared;
- ✓ A 30-minute ion session is identified by **ION-30**, 10 minutes of Zapper by **ZAP-10**, and 3 sessions of 10 minutes of Zapper by **3xZAP-10** in the following grids;
- ✓ Run the Zapper before or after the ionizing bath when they take place on the same day;

- ✓ Here is the summary of the first three months of use of the ionizing bath, including the Zapper:

MONDAY	TUESDAY	WEDNESDAY	THURSDAY	FRIDAY	SATURDAY	SUNDAY
1st MONTH	1 ION-30	2	3	4 ION-30	5	6
7 ION-30	8	9	10 ION-30	11	12	13 ION-30
14	15	16 ION-30	17	18	19 ION-30	20
21	22 ION-30	23	24	25 ION-30	26	27
28 ION-30	29	30	31 ION-30			

MONDAY	TUESDAY	WEDNESDAY	THURSDAY	FRIDAY	SATURDAY	SUNDAY
2nd MONTH	1 ZAP-10 ION-30	2 ZAP-10	3 ZAP-10	4 ZAP-10 ION-30	5 ZAP-10	6 ZAP-10
7 ZAP-10 ION-30	8 ZAP-10	9 ZAP-10	10 ZAP-10 ION-30	11 ZAP-10	12 ZAP-10	13 ZAP-10 ION-30
14 ZAP-10	15 ZAP-10	16 ZAP-10 ION-30	17 ZAP-10	18 ZAP-10	19 ZAP-10 ION-30	20 ZAP-10
21 ZAP-10	22 ZAP-10 ION-30	23 ZAP-10	24 ZAP-10	25 ZAP-10 ION-30	26 ZAP-10	27 ZAP-10
28 ZAP-10 ION-30	29 ZAP-10	30 ZAP-10	31 ZAP-10 ION-30			

MONDAY	TUESDAY	WEDNESDAY	THURSDAY	FRIDAY	SATURDAY	SUNDAY
3rd MONTH	1 3xZAP-10 ION-30	2 3xZAP-10	3 3xZAP-10	4 3xZAP-10 ION-30	5 3xZAP-10	6 3xZAP-10
7 3xZAP-10 ION-30	8 3xZAP-10	9 3xZAP-10	10 3xZAP-10 ION-30	11 3xZAP-10	12 3xZAP-10	13 3xZAP-10 ION-30
14 3xZAP-10	15 3xZAP-10	16 3xZAP-10 ION-30	17 3xZAP-10	18 3xZAP-10	19 3xZAP-10 ION-30	20 3xZAP-10
21 3xZAP-10	22 3xZAP-10 ION-30	23 3xZAP-10	24 3xZAP-10	25 3xZAP-10 ION-30	26 3xZAP-10	27 3xZAP-10
28 3xZAP-10 ION-30	29 3xZAP-10	30 3xZAP-10	31 3xZAP-10 ION-30			

- ✓ After the execution of this initial 3-month protocol, stop the Zapper sessions to let the immune system equilibrate without any intervention;

- ✓ **Resume 1 to 3 sessions of 10 minutes each day** for 1 to 2 weeks to stimulate the immune functions in the event of a resurgence of the following signs:
 - Presence of infection, flu, cold, etc.,
 - Irresistible carbohydrate cravings or daily sugar cravings,
 - Chronic inflammatory pain, excess mucus and secretions,
 - Acidified wounds, itchy skin, and abnormal nervousness.
- ✓ Drink plenty of water and eat little and very healthily while avoiding sugars, carbohydrates, and dairy products to minimize inflammation, metabolic acidity, and the proliferation of harmful parasitic microorganisms in the body;
- ✓ It is possible to combine the sessions with purifying natural products or pure and good quality essential supplements;
- ✓ In case of persistence of certain symptoms or for a particular follow-up, please contact us at info@VitalAxion.com or 1-866-981-0357 to be guided towards other alternative solutions or to be advised on complementary natural methods.